

SLIC Cloning

Reference: Li *et al*, Nature Methods, 2007 “Harnessing homologous recombination in vitro to generate recombinant DNA via SLIC”

Materials:

T4 DNA polymerase (T4 DNAP): NEB (# M0203S)

RecA: NEB (# M0249S)

T4 ligase buffer (containing ATP): NEB (sold together with T4 ligase # M0202S)

100x BSA: NEB (sold together with many enzymes)

Primer design:

Add 30 bps from the vector sequence to 5' of the PCR primers to ensure annealing.

Insert:

PCR amplify the insert in a 50 uL rxn, gel purify product, elute in 30 uL water.

Vector:

Double digest the vector: 2 ug vector in 50 uL rxn, 1-2 hr at desired temperature, gel purify, elute in 30 uL water. Or, PCR amplify the vector backbone – especially if you do not want to rely on restriction sites.

SLIC cloning:

1. T4 DNAP treatment

- Make master T4 RNAP solution for 12 rxns:

60 uL

1.2 uL BSA (10mg/ml from NEB)

12 uL 10x T4 ligase buffer

1 uL T4 DNAP stock

45.8 uL H2O

- 5 uL of master solution + 5 uL eluted DNA
- RT 30min; then stop rxn by adding 1 uL 10 mM dCTP and more to ice

2. Annealing

- Make master RecA solution:
 - 1) Dilute RecA: 1 uL RecA + 5 uL ddH2O
 - 2) Take 1 uL diluted RecA, add 1 uL BSA + 10 uL T4 ligase buffer + 88 uL water
- 6 uL master RecA + 2 uL treated vector + 2 uL treated insert
- 37 C 30 min; then move to ice ready for transformation

3. Transformation

Use 5-10 uL to transform 100 uL Mach1 competent cells. Do not use any short-cut transformation methods.

SLIC Cloning